



Case number CAS-0023761384

Order number WO-33511965 | Visit Date 29/07/2025

Instrument Details

Instrument/Module cobas 6800

Serial Number 1376

Internal Instrument Name

Additional Details

PO

Customer Details

Customer Number 0052027574

Customer Name P'CHOMKLAO HOSPITAL
รพ.พระจอมเกล้า

Street Address 53 ถ.รถไฟ ต.คลองกระเซง

ZIP 76000 - Phetchaburi

Contact Name Supaporn Sae-aeaw

Contact Phone +66864684997

Service Activity Code Remote Support

Purpose of Visit

Related instrument: cobas 6800 sn. 1376.

Location: Phrachomklao hospital.

Issue : Transfer module has changed to error during sample transferring step.

The issue related to the sample pipetter head C.

Possible cause: Sample transfer head C.

Loss of reagent: No wastage reagent.

Action taken: Suggested customer that performed the instrument clean up and then restarted the analyzer. The system was down and could not be run normally. The issue need to be further investigated.

Note: Transfer the case to FSE for fixation and resolution.

Performed Activities

Transfer module changed to error when sample transferring step.

Description

Related instrument: cobas 6800 sn. 1376.

Location: Phrachomklao hospital.

Issue : Transfer module has changed to error during sample transferring step.

The issue related to the sample pipetter head C.

Possible cause: Sample transfer head C.

Loss of reagent: No wastage reagent.

Action taken: Suggested customer that performed the instrument clean up and then restarted the analyzer. The system was down and could not be run normally. The issue need to be further investigated.

Note: Transfer the case to FSE for fixation and resolution.

Time Report

Category	Start	End	Hours	Billing Type
05913616001-Service labour time	29/07/2568 12:30	29/07/2568 14:00	1.5	Free of Charge

Signature or Stamp

Customer
Supaporn Sae-aeaw

Roche Support
Kasama Rakmark

Sign here

By signing this document

- Customer representative confirms they have reviewed the work carried out as described and accepts the system for use.
- Roche Support person confirms the above mentioned work and checks have been completed..

- .

Hotline : 1800-0193-99
Tel : 02 7912200
Line Id : @rcscthailand

[Roche Customer Support Center](#)